


# LeMaitre® Embolectomy Catheters



## ORDERING INFORMATION

Balloon Diameter	Balloon Volume	Diameter	Color	Length	Guidewire	Model #
<b>Single Lumen Embolectomy Catheter</b>						
4.5 mm	0.05 ml	2F	White	40 cm		e1601-24
4.5 mm	0.05 ml	2F	White	60 cm		e1601-26
8.0 mm	0.20 ml	3F	Green	40 cm		e1601-34
8.0 mm	0.20 ml	3F	Green	80 cm		e1601-38
10.5 mm	0.75 ml	4F	Red	40 cm		e1601-44
10.5 mm	0.75 ml	4F	Red	80 cm		e1601-48
13.0 mm	1.50 ml	5F	White	80 cm		e1601-54
13.0 mm	1.50 ml	5F	White	80 cm		e1601-58
13.5 mm	1.60 ml	6F	Blue	80 cm		e1601-68
14.0 mm	1.75 ml	7F	Yellow	80 cm		e1601-78

## Over-the-Wire Embolectomy Catheter

Compatible with standard guidewire sizes. Stopcock and syringe included.

6 mm	0.20 ml	3F	Green	40 cm	.018"	e1651-34
6 mm	0.20 ml	3F	Green	80 cm	.018"	e1651-38
10 mm	0.75 ml	4F	Red	40 cm	.025"	e1651-44
10 mm	0.75 ml	4F	Red	80 cm	.025"	e1651-48
12 mm	1.50 ml	5F Plus	White	40 cm	.035"	e1651-84
12 mm	1.50 ml	5F Plus	White	80 cm	.035"	e1651-88
13 mm	1.60 ml	6F	Blue	40 cm	.035"	e1651-64
13 mm	1.60 ml	6F	Blue	80 cm	.035"	e1651-68
14 mm	1.75 ml	7F	Yellow	80 cm	.038"	e1651-78

## NovaSil™ Silicone Single Lumen Embolectomy Catheter\*

4 mm	0.05 ml	2F	White	60 cm		e1801-26
6 mm	0.20 ml	3F	Green	40 cm		e1801-34
6 mm	0.20 ml	3F	Green	80 cm		e1801-38
9 mm	0.60 ml	4F	Red	40 cm		e1801-44
9 mm	0.60 ml	4F	Red	80 cm		e1801-48
11 mm	1.00 ml	5F	White	40 cm		e1801-54
11 mm	1.00 ml	5F	White	80 cm		e1801-58
13 mm	1.60 ml	6F	Blue	80 cm		e1801-68
14 mm	1.75 ml	7F	Yellow	80 cm		e1801-78

## LeMaitre Catheter Rack

Catheter rack compatible for all single lumen, over-the-wire, and silicone catheters.

Catheter Rack	A3000-00
Wall Mount	A3000-01

These specifications are not intended as a warranty. In the interest of product improvement, these specifications may be changed from time to time without notice. Please consult your sales representative for details.

\*Not made with natural rubber latex.

\*Comparative bench top testing performed on the LeMaitre 4F Single Lumen Latex Embolectomy Catheter and the leading competitor's 4F catheter demonstrated up to 80% stronger balloon bonds and up to 60% stronger balloons.

The following models are not licensed for sale in Canada: all e1801-XX model #'s and e1651-84, e1651-88.



Your Peripheral Vision™

LeMaitre is a registered trademark and NovaSil and Your Peripheral Vision are trademarks of LeMaitre Vascular, Inc.  
©2013 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA. M0005 Rev. L 04/13



### LeMaitre Vascular, Inc.

63 Second Avenue  
Burlington, MA 01803  
USA

T 781 221 2266

F 781 221 2223

### LeMaitre Vascular ULC

5915 Airport Road, Suite 608  
Mississauga, ON L4V 1T1  
Canada

T 905 673 2266

F 905 673 2223

### LeMaitre Vascular GmbH

Otto-Volger-Str. 5a/b  
65843 Sulzbach/Ts.  
Germany

T +49 (0)6196 659230

F +49 (0)6196 527072

### LeMaitre Vascular GK

1F Kubodera Twin Tower Bldg.  
2-9-4 Kudan-minami, Chiyoda-ku  
Tokyo 102-0074  
Japan

T +81 (0)3 5215 5681

F +81 (0)3 5215 5682

### Customer Service

US 800 628 9470

Canada 855 673 2266

www.lemaitre.com

# LeMaitre® Embolectomy Catheters

The World's Strongest Embolectomy Catheter



Your Peripheral Vision™

# LeMaitre<sup>®</sup> Embolectomy Catheters

## STRONGER BOND. STRONGER BALLOON. UNBEATABLE PRICE.\*

### Superior strength.

The significantly stronger LeMaitre Single Lumen Latex Embolectomy Catheter reduces frustration and enhances safety. It features stronger bonds and a stronger balloon at an unbeatable price.

### New low price

How often can you buy a superior product for less money? Not only is the LeMaitre Single Lumen Embolectomy Catheter stronger and less prone to rupture, it costs less per unit than the leading competitor.

### Trusted brand

LeMaitre Embolectomy Catheters have been sold internationally for over 20 years and have been used in over 1 million procedures.

### OVER-THE-WIRE

#### Guidance.

Compatible with standard guidewire to direct catheter to proper position

#### Dye Enhancement.

Inject contrast media to facilitate emboli/thrombi imaging.

#### Irrigation.

Allows for lubrication during catheter insertion and withdrawal

### STRONGER BOND TEST RESULTS\*

(lbf) = pounds of force



**LONGER  
SHELF LIFE**  
Latex: 6 yrs  
Silicone: 5 yrs

### SINGLE LUMEN

60% stronger, rupture-resistant balloons\*  
80% stronger catheter bonds\*

Stronger bonds and balloons allow you to apply more force with less risk.

### SILICONE

The silicone balloon catheter provides a viable option for latex sensitivity.



Your Peripheral Vision™